

# **Product Data Sheet**

## hsa-miR-6090 miRNA Agomir

| Catalog #               | Source      | Reactivity   | Applications  |  |  |  |
|-------------------------|-------------|--|---|--|--|--|
| CMJ1965                 | Synthetic   | н  |   |  |  |  |
| Description             | Synth       | etic miRNA Agomir is u   | sed to regulate the expression of target hsa-miR-6090     |  |  |  |
|                         | mRNA        | ۹.   |   |  |  |  |
| Specificity             | Agom        | ir is chemically-modifie   | d double-strand miRNA mimics which can mimic              |  |  |  |
|                         | matu        | re endogenous miRNAs   | after transfection into cells. They can up-regulate the   |  |  |  |
|                         | endo        | genous miRNA activity  | by utilizing the natural miRNA machinery. Our miRNA       |  |  |  |
|                         | agom        | ir is designed to mimic  | mature miRNAs, and chemically modified to increase        |  |  |  |
|                         | their       | stability and activity.  |   |  |  |  |
| Form                    | Lyoph       | ilized powder  |   |  |  |  |
| Gene Symbol             | hsa-m       | niR-6090   |   |  |  |  |
| Accession No.           | MIMA        | AT0023715  |   |  |  |  |
| Components              | This s      | This synthetic miRNA is based on the mature miRNA sequence. The antisense strand |   |  |  |  |
|                         | of the      | e agomir has 2 phospho   | rothioates at the 5' end, 4 phosphorothioates, 1          |  |  |  |
|                         | chole       | sterol group at the 3' e   | nd, and full-length nucleotide 2'-methoxy modification.   |  |  |  |
|                         | lt exh      | ibits enhanced cellular  | uptake, stability and regulatory activity in vivo.        |  |  |  |
| Directions for <b>l</b> | Jse Briefly | y centrifuge tubes cont  | aining miRNA agomir to ensure that the miRNA pellet is    |  |  |  |
|                         | locate      | ed at the bottom of the  | tube. Dissolve miRNA agomir to a convenient stock         |  |  |  |
|                         | conce       | entration using the reco   | mmended volume of DEPC H2O (or RNase-free water).         |  |  |  |
|                         | For ex      | ample: dissolve 5 nmo  | miRNA agomir to 20 $\mu M$ using 250 $\mu l$ DEPC H2O (or |  |  |  |
|                         | RNase       | e-free water). Pipette th  | e solution up and down 3-5 times (or vortex briefly).     |  |  |  |
|                         | Briefly     | y centrifuge tubes cont  | aining miRNA agomir to ensure that the solution is        |  |  |  |
|                         | collec      | ted at the bottom of th  | e tube. Aliquot the miRNA agomir into small volumes       |  |  |  |
| Application kow         |             | loctorn blot IH Immunobi   | stochomistry IE Immunofluorosconco EC Elow sytomatry IC   |  |  |  |

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC-Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb-Rabbit, S- Sheep, Z- Zebrafish

### **COHESION BIOSCIENCES LIMITED**

| WEB                 | ORDER                 | SUPPORT                     | CUSTOM                 |
|---------------------|-----------------------|-----------------------------|------------------------|
| www.cohesionbio.com | order@cohesionbio.com | techsupport@cohesionbio.com | custom@cohesionbio.com |



## **Product Data Sheet**

and store at  $\leq$  -20°C. miRNA agomir is stable (for 6 months under the specified storage condition). For best results, use in 3 months and limit freeze-thaw events for each tube no more than five times.

**Storage/Stability** Shipped at 4 °C. Store at -20 °C for one year. Avoid freeze-thaw cycles after reconstitution.

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC-Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb-Rabbit, S- Sheep, Z- Zebrafish

#### **COHESION BIOSCIENCES LIMITED**

 WEB
 ORDER
 SUPPORT
 CUSTOM

 www.cohesionbio.com
 order@cohesionbio.com
 techsupport@cohesionbio.com
 custom@cohesionbio.com