

Anti-TAU Antibody (Detection)

| Catalog # | Source | Reactivity | Applications |
|---------------------------|---|------------|--------------|
| CAP1118-D | Mouse | Human | SE |
| Description | Mouse monoclonal detection antibody to TAU | | |
| Immunogen | Recombinant protein of human TAU | | |
| Purification | The antibody was purified by Protein A/G affinity chromatography. | | |
| Specificity | Recognizes native and recombinant human TAU protein | | |
| Clonality | Monoclonal | | |
| Conjugation | | | |
| Form | Mouse IgG2b. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide. | | |
| Dilution | | | |
| Gene Symbol | MAPT | | |
| Alternative Names | MAPTL; MTBT1; TAU; Microtubule-associated protein tau; Neurofibrillary tangle protein; Paired helical filament-tau; PHF-tau | | |
| Entrez Gene | 4137 | | |
| SwissProt | P10636-8 | | |
| Directions for Use | Recommended pair for sandwich immunoassay - Capture: CAP1118-C | | |
| Storage/Stability | Shipped at 4° C. Upon delivery aliquot and store at -20° C for 1 year. Avoid freeze/thaw cycles. | | |

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com