

Tyr-CRF (Ovine)

Catalog #	Source	Reactivity	Applications
CCP3003	Synthetic		
Description	Peptide to Tyr-CRF (Ovine)		
Form	Lyophilized powder		
CAS Number	83930-34-1		
Molecular Formula	C214H348N60O65S		
Molecular Weight	4833.59		
Purity	> 95%		
Chemical Structure	Tyr - Ser - Gln - Glu - Pro - Pro - Ile - Ser - Leu - Asp - Leu - Thr - Phe - His - Leu - Leu - Arg - Glu - Val - Leu - Glu - Met - Thr - Lys - Ala - Asp - Gln - Leu - Ala - Gln - Gln - Ala - His - Ser - Asn - Arg - Lys - Leu - Leu - Asp - Ile - Ala - NH2		
Storage/Stability Structure	Shipped at 4°C. Store at -20°C for one year.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com