

AGRP (87-132) (Human)

Catalog #	Source	Reactivity	Applications
CCP3017	Synthetic		
Description	Peptide to AGRP (87-132) (Human)		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C ₂₁₉ H ₃₃₉ N ₆₅ O ₆₃ S ₁₁		
Molecular Weight	5243.17		
Purity	> 95%		
Chemical Structure	Ac - Cys - Val - Arg - Leu - His - Glu - Ser - Cys - Leu - Gly - Gln - Gln - Val - Pro - Cys - Cys - Asp - Pro - Cys - Ala - Thr - Cys - Tyr - Cys - Arg - Phe - Phe - Asn - Ala - Phe - Cys - Tyr - Cys - Arg - Lys - Leu - Gly - Thr - Ala - Met - Asn - Pro		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com