

**Adrenomedullin (1-50) (Rat)**

Catalog #	Source	Reactivity	Applications
CCP3026	Synthetic		
<b>Description</b>	Peptide to Adrenomedullin (1-50) (Rat)		
<b>Form</b>	Lyophilized powder		
<b>CAS Number</b>			
<b>Molecular Formula</b>	C248H381N77O75S5		
<b>Molecular Weight</b>	5729.50		
<b>Purity</b>	> 95%		
<b>Chemical Structure</b>	Tyr - Arg - Gln - Ser - Met - Asn - Gln - Gly - Ser - Arg - Ser - Thr - Gly - Cys - Arg - Phe - Gly - Thr - Cys - Thr - Met - Gln - Lys - Leu - Ala - His - Gln - Ile - Tyr - Gln - Phe - Thr - Asp - Lys - Asp - Lys - Asp - Gly - Met - Ala - Pro - Arg - Asn		
<b>Storage/Stability</b>	Shipped at 4°C. Store at -20°C for one year.		
<b>Structure</b>			

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

**COHESION BIOSCIENCES LIMITED**

**WEB**  
[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)