

rec Acetyl-Eglin c

Catalog #	Source	Reactivity	Applications
CCP3030	Synthetic		
Description	Peptide to rec Acetyl-Eglin c		
Form	Lyophilized powder		
CAS Number	96380-69-7		
Molecular Formula	C375H552N96O108		
Molecular Weight	8133.16		
Purity	> 95%		
Chemical Structure	Ac - Thr - Glu - Phe - Gly - Ser - Glu - Leu - Lys - Ser - Phe - Pro - Glu - Val - Val - Gly - Lys - Thr - Val - Asp - Gln - Ala - Arg - Glu - Tyr - Phe - Thr - Leu - His - Tyr - Pro - Gln - Tyr - Asp - Val - Tyr - Phe - Leu - Pro - Glu - Gly - Ser - Pro		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com