

NTproBNP (1-76)

Catalog #	Source	Reactivity	Applications
CCP3032	Synthetic		
Description	Peptide to NTproBNP (1-76)		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C364H596N114O116S		
Molecular Weight	8457.57		
Purity	> 95%		
Chemical Structure	His - Pro - Leu - Gly - Ser - Pro - Gly - Ser - Ala - Ser - Asp - Leu - Glu - Thr - Ser - Gly - Leu - Gln - Glu - Gln - Arg - Asn - His - Leu - Gln - Gly - Lys - Leu - Ser - Glu - Leu - Gln - Val - Glu - Gln - Thr - Ser - Leu - Glu - Pro - Leu - Gln - Glu		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com