

hsa-miR-1255b-2-3p miRNA Inhibitor

Catalog #	Source	Reactivity	Applications
CIH1893	Synthetic	H	
Description	Synthetic miRNA inhibitors are used to inhibit hsa-miR-1255b-2-3p by transfection.		
Specificity	miRNA Inhibitors are chemically synthesized, complementary, antisense single-stranded oligonucleotides to their target, endogenous mature miRNA. It effectively prevents the target miRNA to bind to normal cellular binding sites. Transfection of inhibitors followed by downstream gene expression analysis or phenotypic analysis, is performed to elucidate the targets and roles of particular miRNAs.		
Form	Lyophilized powder		
Gene Symbol	hsa-miR-1255b-2-3p		
Accession No.	MIMAT0022725		
Components	This synthetic miRNA is based on the mature miRNA sequence. It does not contain the full precursor miRNA stem-loop.		
Directions for Use	We recommend re-suspending the lyophilized synthetic miRNA using DNase and RNase-free ddH ₂ O. To make a 100 uM stock, dissolve the lyophilized powder using 100 ul of ddH ₂ O.		
Storage/Stability	Shipped at 4 °C. Store at -20 °C for one year. Avoid freeze-thaw cycles after reconstitution.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com