

Product Data Sheet

mmu-miR-344d-2-5p miRNA Inhibitor

Catalog #	Source	Reactivity	Applications
CIM0822	Synthetic	M	
Description	Synthetic miRNA inhibitors are used to inhibit mmu-miR-344d-2-5p by transfection.		
Specificity	miRNA Inhibitors are chemically synthesized, complementary, antisense single-stranded oligonucleotides to their target, endogenous mature miRNA. It effectively prevents the target miRNA to bind to normal cellular binding sites. Transfection of inhibitors followed by downstream gene expression analysis or phenotypic analysis, is performed to elucidate the targets and roles of particular miRNAs.		
Form	Lyophilized powder		
Gene Symbol	mmu-miR-344d-2-5p		
Accession No.	MIMAT0014961		
Components	This synthetic miRNA is based on the mature miRNA sequence. It does not contain the full precursor miRNA stem-loop.		
Directions for Use	We recommend re-suspending the lyophilized synthetic miRNA using DNase and RNase-free ddH ₂ O. To make a 100 uM stock, dissolve the lyophilized powder using 100 ul of ddH ₂ O.		
Storage/Stability	Shipped at 4 °C. Store at -20 °C for one year. Avoid freeze-thaw cycles after reconstitution.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com