

ELISA Coating Buffer

Catalog #	Source	Reactivity	Applications
CRG1111			ELISA
Description	ELISA Coating Buffer is intended for use as a plate coating buffer for developing a sandwich ELISA.		
Form	Liquid		
Directions for Use	This product is provided as a 1X buffer solution that is ready for use. No dilution is necessary.		
Components	The buffer is a 50 mM carbonate buffer, pH 9.4.		
Storage/Stability	Store at 4°C for one year.		

COHESION BIOSCIENCES LIMITED