

## ELISA Stop Solution (For TMB)

Catalog #	Source	Reactivity	Applications
CRG1112			ELISA
<b>Description</b>	ELISA Stop Solution for TMB Substrate is a ready-to-use solution for colorimetric HRP-based ELISA detection.		
<b>Form</b>	Liquid		
<b>Directions for Use</b>	<p>(a) After incubation with horseradish peroxidase (HRP) conjugated secondary antibody, wash plates with PBS/0.05% Tween-20.</p> <p>(b) Add 100 <math>\mu</math>l of TMB solution/well. Incubate for a specific amount of time or until desired signal is achieved.</p> <p>(c) Add 100 <math>\mu</math>l of ELISA Stop Solution.</p> <p>(d) Read at Absorbance 450 nm.</p>		
<b>Storage/Stability</b>	Store at 4°C for one year.		

### COHESION BIOSCIENCES LIMITED