

Recombinant Human LAIR1 Protein

Catalog #	Source	Reactivity	Applications
CRP1872	Human cells	Human	E, WB, SDS-PAGE, MS

Description Recombinant Human LAIR1 Protein is produced by our Mammalian expression system and the target gene encoding Gln22-His163 is expressed with a 6His tag at the C-terminus.

Form Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Gene Symbol LAIR1

Alternative Names CD305; Leukocyte-associated immunoglobulin-like receptor 1; LAIR-1; hLAIR1; CD305

Entrez Gene 3903 (Human)

SwissProt Q6GTX8 (Human)

Purity Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Chemical Structure QEEDLPRPSI SAEPGTVIPL GSHVTFVCRG PVGVQTFRLE RESRSTYNDT EDVSQASPSE SEARFRIDSV SEGNAGPYRC IYYKPPKWSE QSDYLELLVK ETSGGPDSPD TEPGSSAGPT QRPSDNSHNE HAPASQGLKA EHV DHHHHHHH

Quality Control Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Directions for Use Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Storage/Stability Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com

Product Data Sheet

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com