

Recombinant Human LALBA Protein

Catalog #	Source	Reactivity	Applications
CRP1874	Human cells	Human	E, WB, SDS-PAGE, MS
Description	Recombinant Human LALBA Protein is produced by our mammalian expression system and the target gene encoding Lys20-Leu142 is expressed with a 6His tag at the C-terminus.		
Form	Liquid in a 0.2 μM filtered solution of 20mM Tris HCl, 150mM NaCl, pH8.0.		
Gene Symbol	LALBA		
Alternative Names	LYZL7; Alpha-lactalbumin; Lactose synthase B protein; Lysozyme-like protein 7		
Entrez Gene	3906 (Human)		
SwissProt	P00709 (Human)		
Purity	Greater than 95% as determined by reducing SDS-PAGE.		
Chemical Structure	KQFTKCELSQ LLKDIDGYGG IALPELICTM FHTSGYDTQA IVENNESTEY GLFQISNKLW CKSSQVPQSR NICDISCDKF LDDDITDDIM CAKKILDIKG IDYWLAHKAL CTEKLEQWLC EKLVDHHHHH H		
Quality Control	Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.		
Directions for Use	Always centrifuge tubes before opening. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com