

Recombinant Human UBE2I Protein

Catalog #	Source	Reactivity	Applications
CRP2482	E. coli	Human	E, WB, SDS-PAGE, MS
Description	Recombinant Human UBE2I Protein is produced by our E. coli expression system and the target gene encoding Met1-Ser158 is expressed with a GST tag at the N-terminus.		
Form	Liquid in a 0.2 μM filtered solution of 50mM HEPES, 150mM NaCl, pH 7.5.		
Gene Symbol	UBE2I		
Alternative Names	UBC9; UBCE9; SUMO-conjugating enzyme UBC9; SUMO-protein ligase; Ubiquitin carrier protein 9; Ubiquitin carrier protein I; Ubiquitin-conjugating enzyme E2 I; Ubiquitin-protein ligase I; p18		
Entrez Gene	7329 (Human)		
SwissProt	P63279 (Human)		
Purity	Greater than 95% as determined by reducing SDS-PAGE.		
Chemical Structure	<p>MSPILGYWKI KGLVQPTRL LLEYLEEKYEE HLYERDEGDK WRNKKFELGL EFPNLPYYID</p> <p>GDVKLTQSMA IIRYIADKHN MLGGCPKERA EISMLEGAVL DIRYGVSRIA YSKDFETLKV</p> <p>DFLSKLP EML KMFEDRLCHK TYLNGDHVTH PDFMLYDALD VVLYMDPMCL DAFPKLVCFK</p> <p>KRIEAIPIQID KYLKSSKYIA WPLQGWQATF GGGDHPPKSD LVPRGSHMSG IALSRLAQER</p> <p>KAWRKDHPFG FVAVPTKNPD GTMNL MNWEC AIPGKKGTPW EGGLFKLRML</p> <p>FKDDYPSSPP KCKFEPPLFH PNVYPSGTVC LSILEEDKDW RPAITIKQIL LGIQELLNEP</p> <p>NIQDPAQAEA YTIYCQNRVE YEKRVRAQAK KFAPS</p>		
Quality Control	Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.		
Directions for Use	Always centrifuge tubes before opening. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com

Product Data Sheet

Storage/Stability Store it at -20°C to -80°C for one year.

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com