

Recombinant Human UBE2M Protein

| Catalog # | Source | Reactivity | Applications |
|---------------------------|--|------------|---------------------|
| CRP2487 | E. coli | Human | E, WB, SDS-PAGE, MS |
| Description | Recombinant Human UBE2M Protein is produced by our E. coli expression system and the target gene encoding Met1-Lys183 is expressed. | | |
| Form | Liquid in a 0.2 μM filtered solution of 50mM HEPES, 2mM DTT, 150mM NaCl, 10% Glycerol, pH 7.5. | | |
| Gene Symbol | UBE2M | | |
| Alternative Names | UBC12; NEDD8-conjugating enzyme Ubc12; NEDD8 carrier protein; NEDD8 protein ligase; Ubiquitin-conjugating enzyme E2 M | | |
| Entrez Gene | 9040 (Human) | | |
| SwissProt | P61081 (Human) | | |
| Purity | Greater than 95% as determined by reducing SDS-PAGE. | | |
| Chemical Structure | MIKLFSLKQQ KKEESAGGT KGSSKKASAA QLRIQKDINE LNLPKTCDIS FSDPDDLNF KLVICPDEGF YKSGKFVFSF KVGQGYPHDP PKVKCETMVY HPNIDLEGNV CLNILREDWK PVLINSIIY GLQYLFLEPN PEDPLNKEAA EVLQNNRRLF EQNVQRSMRG GYIGSTYFER CLK | | |
| Quality Control | Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test. | | |
| Directions for Use | Always centrifuge tubes before opening. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. | | |
| Storage/Stability | Store it at -20°C to -80°C for one year. | | |

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com