

## Product Data Sheet

## Mouse Anti-Human IgE - HRP

| Catalog #                | Source                                                                                                                                                   | Reactivity | Applications  |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| CSA2184                  | Mouse                                                                                                                                                    | H          | E, WB, IH, IC |
| <b>Description</b>       | Mouse Monoclonal Secondary Antibody to Human IgE HRP labeled                                                                                             |            |               |
| <b>Immunogen</b>         | Human IgE                                                                                                                                                |            |               |
| <b>Purification</b>      | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.                                         |            |               |
| <b>Specificity</b>       | By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgE. No antibody was detected against non immunoglobulin serum proteins. |            |               |
| <b>Clonality</b>         | Monoclonal                                                                                                                                               |            |               |
| <b>Conjugation</b>       | HRP                                                                                                                                                      |            |               |
| <b>Form</b>              | 1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Proclin                                               |            |               |
| <b>Dilution</b>          | E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)                                                                    |            |               |
| <b>Storage/Stability</b> | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                                         |            |               |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

WEB  
[www.cohesionbio.com](http://www.cohesionbio.com)

ORDER  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

SUPPORT  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

CUSTOM  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)