

Hexokinase

| Catalog # | Source | Reactivity | Applications |
|--------------------------|--|------------|--------------|
| CZM1011 | Microorganism | | |
| Description | Hexokinase from Microorganism | | |
| Specificity | >400 U/mg protein, >100 U/mg powder | | |
| Form | White Lyophilized powder | | |
| Alternative Names | EC 2.7.1.1; HK; hexokinase type IV glucokinase; hexokinase D; hexokinase type IV; hexokinase (phosphorylating); ATP-dependent hexokinase; glucose ATP phosphotransferase | | |
| CAS Number | 9001-51-8 | | |
| Molecular Weight | | | |
| Purity | >90% | | |
| Storage/Stability | Shipped at 4° C. Store at -20° C for 2 years. | | |

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com