

## Citrate Lyase

Catalog #	Source	Reactivity	Applications
CZM1059	Microorganism		
<b>Description</b>	Citrate Lyase from Microorganism		
<b>Specificity</b>	>0.25 U/mg		
<b>Form</b>	Slightly beige lyophilizate		
<b>Alternative Names</b>	EC 4.1.3.6; hACL2		
<b>CAS Number</b>	9012-83-3		
<b>Molecular Weight</b>			
<b>Purity</b>	>90%		
<b>Storage/Stability</b>	Shipped at 4° C. Store at -20° C for 2 years.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)