

Citrate Synthase

Catalog #	Source	Reactivity	Applications
CZM1060	Porcine heart		
Description	Citrate Synthase from Porcine heart		
Specificity	>100 U/mg protein		
Form	Lyophilized powder		
Alternative Names	EC 4.1.3.7; CS; Citrate (Si)-synthase; (R)-citric synthase; Citrate oxaloacetate-lyase [[pro-3S)-CH ₂ COO-->acetyl-CoA]		
CAS Number	9027-96-7		
Molecular Weight			
Purity	>90%		
Storage/Stability	Shipped at 4° C. Store at -20° C for 2 years.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com