

Hyaluronidase

Catalog #	Source	Reactivity	Applications
CZM1069	Bovine Testis		

Description	Hyaluronidase from Bovine Testis
Specificity	>350 U/mg powder
Form	Yellowish white powder
Alternative Names	EC 3.2.1.35; hyaluronoglucosaminidase; hyaluronidase; hyaluronoglucosidase; chondroitinase; chondroitinase I; hyaluronate 4-glycanohydrolase
CAS Number	37326-33-3
Molecular Weight	
Purity	>90%
Storage/Stability	Shipped at 4° C. Store at -20° C for 2 years.

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com