

Alanine Aminotransferase

Catalog #	Source	Reactivity	Applications
CZM1072	Porcine heart		
Description	Alanine Aminotransferase from Porcine heart		
Specificity	>1 U/mg protein, >1 U/mg powder		
Form	White lyophilized powder		
Alternative Names	EC 2.6.1.2; ALT; ALT; alanine aminotransferase; ALAT; glutamic-pyruvic transaminase; glutamic-alanine transaminase; GPT; β -alanine aminotransferase; alanine- α -ketoglutarate aminotransferase; alanine-pyruvate aminotransferase; glutamic acid-pyruvic acid transaminase; glutamic-pyruvic aminotransferase; L-alanine aminotransferase; L-alanine transaminase; L-alanine- α -ketoglutarate aminotransferase; pyruvate transaminase; pyruvate-alanine aminotransferase; pyruvate-glutamate transaminase		
CAS Number	9000-86-6		
Molecular Weight			
Purity	>90%		
Storage/Stability	Shipped at 4° C. Store at -20° C for 2 years.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com