

**RSAD1 Blocking Peptide**

| Catalog #                 | Source  | Reactivity | Applications |
|---------------------------|---|------------|--------------|
| CBP5022                   | Synthetic   | H, M, B    | BL           |
| <b>Description</b>        | The peptide is used to block Anti-RSAD1 Antibody (#CPA5022) reactivity.   |            |              |
| <b>Form</b>               | Lyophilized powder  |            |              |
| <b>Gene Symbol</b>        | RSAD1   |            |              |
| <b>Alternative Names</b>  | Radical S-adenosyl methionine domain-containing protein 1, mitochondrial;<br>Oxygen-independent coproporphyrinogen-III oxidase-like protein RSAD1                                 |            |              |
| <b>Entrez Gene</b>        | 55316 (Human); 237926 (Mouse)   |            |              |
| <b>SwissProt</b>          | Q9HA92 (Human); Q5SUV1 (Mouse)  |            |              |
| <b>Purity</b>             | >85%  |            |              |
| <b>Quality Control</b>    | The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.   |            |              |
| <b>Directions for Use</b> | Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |            |              |
| <b>Storage/Stability</b>  | Shipped at 4°C. Store at -20°C for one year.  |            |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

**COHESION BIOSCIENCES LIMITED**

**WEB**  
[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)