

MLL1 Blocking Peptide

Catalog #	Source	Reactivity	Applications
CBP6103	Synthetic	H, M, R, B, D, Mk, P	BL
Description	The peptide is used to block Anti-MLL1 Antibody (#CPA6103) reactivity.		
Form	Lyophilized powder		
Gene Symbol	MLL		
Alternative Names	ALL1; CXXC7; HRX; HTRX; MLL; MLL1; TRX1; Histone-lysine N-methyltransferase 2A; Lysine N-methyltransferase 2A; ALL-1; CXXC-type zinc finger protein 7; Myeloid/lymphoid or mixed-lineage leukemia; Myeloid/lymphoid or mixed-lineage leukemia protein 1; Tritho		
Entrez Gene	4297 (Human); 214162 (Mouse)		
SwissProt	Q03164 (Human); P55200 (Mouse)		
Purity	>85%		
Quality Control	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.		
Directions for Use	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com