

mmu-miR-130b-5p miRNA Agomir

Catalog #	Source	Reactivity	Applications
CMN0406	Synthetic	M	
Description	Synthetic miRNA Agomir is used to regulate the expression of target mmu-miR-130b-5p mRNA.		
Specificity	Agomir is chemically-modified double-strand miRNA mimics which can mimic mature endogenous miRNAs after transfection into cells. They can up-regulate the endogenous miRNA activity by utilizing the natural miRNA machinery. Our miRNA agomir is designed to mimic mature miRNAs, and chemically modified to increase their stability and activity.		
Form	Lyophilized powder		
Gene Symbol	mmu-miR-130b-5p		
Accession No.	MIMAT0004583		
Components	This synthetic miRNA is based on the mature miRNA sequence. The antisense strand of the agomir has 2 phosphorothioates at the 5' end, 4 phosphorothioates, 1 cholesterol group at the 3' end, and full-length nucleotide 2'-methoxy modification. It exhibits enhanced cellular uptake, stability and regulatory activity in vivo.		
Directions for Use	Briefly centrifuge tubes containing miRNA agomir to ensure that the miRNA pellet is located at the bottom of the tube. Dissolve miRNA agomir to a convenient stock concentration using the recommended volume of DEPC H ₂ O (or RNase-free water). For example: dissolve 5 nmol miRNA agomir to 20 μM using 250 μl DEPC H ₂ O (or RNase-free water). Pipette the solution up and down 3-5 times (or vortex briefly). Briefly centrifuge tubes containing miRNA agomir to ensure that the solution is collected at the bottom of the tube. Aliquot the miRNA agomir into small volumes		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com

Product Data Sheet

and store at $\leq -20^{\circ}\text{C}$. miRNA agomir is stable (for 6 months under the specified storage condition). For best results, use in 3 months and limit freeze-thaw events for each tube no more than five times.

Storage/Stability

Shipped at 4°C . Store at -20°C for one year. Avoid freeze-thaw cycles after reconstitution.

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com