

Recombinant Human CDK2 Protein

Catalog #	Source	Reactivity	Applications
CRP1269	E. coli	Human	E, WB, SDS-PAGE, MS
Description	Recombinant Human CDK2 Protein is produced by our E.coli expression system and the target gene encoding Met1-Leu298 is expressed with a 6His tag at the N-terminus.		
Form	Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 200mM NaCl, 1mM DTT, 40% Glycerol, pH 8.0.		
Gene Symbol	CDK2		
Alternative Names	CDKN2; Cyclin-dependent kinase 2; Cell division protein kinase 2; p33 protein kinase		
Entrez Gene	1017 (Human)		
SwissProt	P24941 (Human)		
Purity	Greater than 95% as determined by reducing SDS-PAGE.		
Chemical Structure	<p>MGSSHHHHHH SGLVPRGSH MENFQKVEKI GEGTYGVVYK ARNKLTGEVV ALKKIRLDTE</p> <p>TEGVPSTAIR EISLLKELNH PNIVKLLDVI HTENKLYLVF EFLHQDLKKF MDASALTGIP</p> <p>LPLIKSYLFQ LLQGLAFCHS HRVLHRDLKP QNLLINTEGA IKLADFGGLAR AFGVPVRTYT</p> <p>HEVVTLWYRA PEILLGCKYY STAVDIWSLG CIFAEMVTRR ALFPGDSEID QLFRIFRITLG</p> <p>TPDEVVWPGV TSMPDYKPSF PKWARQDFSK VVPLDEDGR SLLSQMLHYD PNKRISAKAA</p> <p>LAHPFFQDVT KPVPHLRL</p>		
Quality Control	Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.		
Directions for Use	Always centrifuge tubes before opening. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com