

Recombinant Human CHRN3 Protein

| Catalog # | Source | Reactivity | Applications |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------------|
| CRP1303 | Human cells | Human | E, WB, SDS-PAGE, MS |
| Description | Recombinant Human CHRN3 Protein is produced by our mammalian expression system and the target gene encoding Ile25-Leu232 is expressed with a 6His tag at the C-terminus. | | |
| Form | Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4. | | |
| Gene Symbol | CHRN3 | | |
| Alternative Names | Neuronal acetylcholine receptor subunit beta-3 | | |
| Entrez Gene | 1142 (Human) | | |
| SwissProt | Q05901 (Human) | | |
| Purity | Greater than 95% as determined by reducing SDS-PAGE. | | |
| Chemical Structure | IAENEDALLR HLFQGYQKWV RPVLHSNDTI KVFYGLKISQ LVDVDEKNQL MTTNVWLKQE WTDHKLRWNP DDYGGIHSIK VPSESLWLDP IVLFENADGR FEGSLMTKVI VKSNGTVVWT PPASYKSSCT MDVTFFPFDR QNCSMKFGSW TYDGTMVLDI LINENVDRKD FFDNGEWEL NAKGMKGNRR DGVYSYPFIT YSFVLRRLLD HHHHHH | | |
| Quality Control | Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test. | | |
| Directions for Use | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. | | |
| Storage/Stability | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. | | |

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com

Product Data Sheet

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com