

Recombinant Mouse CXCL15 Protein

| Catalog # | Source | Reactivity | Applications |
|---------------------------|--|------------|---------------------|
| CRP1417 | Human cells | Mouse | E, WB, SDS-PAGE, MS |
| Description | Recombinant Mouse CXCL15 Protein is produced by our mammalian expression system and the target gene encoding Gln26-Ala167 is expressed with a 6His tag at the C-terminus. | | |
| Form | Lyophilized from a 0.2 μM filtered solution of 20mM Tris, 150mM NaCl, pH 8.0. | | |
| Gene Symbol | CXCL15 | | |
| Alternative Names | Scyb15; C-X-C motif chemokine 15; Lungkine; Small-inducible cytokine B15 | | |
| Entrez Gene | 20309 (Mouse) | | |
| SwissProt | Q9WVL7 (Mouse) | | |
| Purity | Greater than 95% as determined by reducing SDS-PAGE. | | |
| Chemical Structure | <p>QELRCLCIQE HSEFIPLKLI KNIMVIFETI YCNRKEVIAV PKNQSMICLD PDAPWVKATV</p> <p>GPITNRFLE DLKQKEFPPA MKLLYSVEHE KPLYLSFGRP ENKRIFPFPI RETSRHFADL</p> <p>AHNSDRNFLR DSSEVSLTGS DAVDHHHHHH</p> | | |
| Quality Control | Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test. | | |
| Directions for Use | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. | | |
| Storage/Stability | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. | | |

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com