

## Recombinant Human KIR2DL4 Protein

Catalog #	Source	Reactivity	Applications
CRP1847	Human cells	Human	E, WB, SDS-PAGE, MS
<b>Description</b>	Recombinant Human KIR2DL4 Protein is produced by our mammalian expression system and the target gene encoding Trp22-His242 is expressed with a 6His tag at the C-terminus.		
<b>Form</b>	Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.		
<b>Gene Symbol</b>	KIR2DL4		
<b>Alternative Names</b>	CD158D; KIR103AS; Killer cell immunoglobulin-like receptor 2DL4; CD158 antigen-like family member D; G9P; Killer cell inhibitory receptor 103AS; KIR-103AS; MHC class I NK cell receptor KIR103AS; CD158d		
<b>Entrez Gene</b>	3805 (Human)		
<b>SwissProt</b>	Q99706 (Human)		
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE.		
<b>Chemical Structure</b>	WAHVGGQDKP FCSAWPSAVV PQGGHVTLRC HYRRGFNIPT LYKKDGVVPV ELYNRIFWNS FLISPVTPAH AGTYRCRGFH PHSPTWSAP SNPLVIMVTG LYEKPSLTAR PGPTVRTGEN VTLSCSSQSS FDIYHLSREG EAHELRLPAV PSINGTFQAD FPLGPATHGE TYRCFGSFHG SPYEWSDASD PLPVSVTGNP SSSWSPSTEP SFKTGIARHL HVDHHHHHHH		
<b>Quality Control</b>	Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.		
<b>Directions for Use</b>	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
<b>Storage/Stability</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature		
<b>Application key:</b>	E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference		
<b>Species reactivity key:</b>	H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish		

### COHESION BIOSCIENCES LIMITED

# Product Data Sheet

for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days.

Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

**WEB**

[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**

[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)