

Recombinant Human LIPF Protein

Catalog #	Source	Reactivity	Applications
CRP1906	Human cells	Human	E, WB, SDS-PAGE, MS
Description	Recombinant Human LIPF Protein is produced by our Mammalian expression system and the target gene encoding Leu20-Lys398 is expressed with a 6His tag at the C-terminus.		
Form	Liquid in a 0.2 μM filtered solution of 25mM Tris-HCl, 100mM glycine, 10% Glycerol, pH 7.3.		
Gene Symbol	LIPF		
Alternative Names	Gastric triacylglycerol lipase; GL; Gastric lipase		
Entrez Gene	8513 (Human)		
SwissProt	P07098 (Human)		
Purity	Greater than 95% as determined by reducing SDS-PAGE.		
Chemical Structure	<p>LFGKLHPGSP EVTMNISQMI TYWGYPNEEY EVVTEDGYIL EVNRIPYGKK NSGNTGQRPV</p> <p>VFLQHGLLAS ATNWISNLPN NSLAFILADA GYDVWLGNSR GNTWARRNLY YSPDSVEFWA</p> <p>FSFDEMAKYD LPATIDFIVK KAGQKQLHYV GHSQGTIGF IAFSTNPSLA KRIKTFYALA</p> <p>PVATVKYTKS LINKLRFVPQ SLFKFIFGDK IFYPHNFFDQ FLATEVCSRE MLNLLCSNAL</p> <p>FIICGFDSKN FNTSRLDVYL SHNPAGTSVQ NMFHWTQAVK SGKFAQYDWG</p> <p>SPVQNRMHYD QSQPPYYNVT AMNVPIAVWN GGDLLADPQ DVGLLLPKLP NLIYHKEIPF</p> <p>YNHLDFIWAM DAPQEVYNDI VSMISEDKKV DHHHHHH</p>		
Quality Control	Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.		
Directions for Use	Always centrifuge tubes before opening. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com

Product Data Sheet

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com