

## Recombinant Human NPTN Protein

Catalog #	Source	Reactivity	Applications
CRP2021	Human cells	Human	E, WB, SDS-PAGE, MS
<b>Description</b>	Recombinant Human NPTN Protein is produced by our Mammalian expression system and the target gene encoding Gln29-His220 is expressed with a 6His tag at the C-terminus.		
<b>Form</b>	Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.		
<b>Gene Symbol</b>	NPTN		
<b>Alternative Names</b>	SDFR1; SDR1; Neuroplastin; Stromal cell-derived receptor 1; SDR-1		
<b>Entrez Gene</b>	27020 (Human)		
<b>SwissProt</b>	Q9Y639 (Human)		
<b>Purity</b>	Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.		
<b>Chemical Structure</b>	<p>QNEPRIVTSE EVIIRDSPVL PVTLQCNLTS SSHTLTYSYW TKNGVELSAT RKNASNMEYR</p> <p>INKPRAEDSG EYHCVYHFVS APKANATIEV KAAPDITGHK RSENKNEGQD ATMYCKSVGV</p> <p>PHPDWIWRKK ENGMPMDIVN TSGRFFIINK ENYTELNIVN LQITEDPGEY ECNATNAIGS</p> <p>ASVVTVLRVR SHVDHHHHHH</p>		
<b>Quality Control</b>	Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.		
<b>Directions for Use</b>	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
<b>Storage/Stability</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com

# Product Data Sheet

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

**WEB**

[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**

[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)