

Recombinant Human PBLD Protein

Catalog #	Source	Reactivity	Applications
CRP2057	E. coli	Human	E, WB, SDS-PAGE, MS
Description	Recombinant Human PBLD Protein is produced by our E. coli expression system and the target gene encoding Met1-Ala288 is expressed with a 6His tag at the N-terminus.		
Form	Liquid in a 0.2 µM filtered solution of 20mM Tris, 100mM NaCl, 1mM DTT, 30% Glycerol, pH 8.0.		
Gene Symbol	PBLD		
Alternative Names	MAWBP; Phenazine biosynthesis-like domain-containing protein; MAWD-binding protein; MAWDBP; Unknown protein 32 from 2D-page of liver tissue		
Entrez Gene	64081 (Human)		
SwissProt	P30039 (Human)		
Purity	Greater than 95% as determined by reducing SDS-PAGE.		
Chemical Structure	<p>MGSSHHHHHH SSGLVPRGSH MKLPIFIADA FTARAFRGNP AAVCLLENEL DEDMHQKIAR</p> <p>EMNLSETAFI RKLHPTDNFA QSSCFGLRWF TPASEVPLCG HATLASAAVL FHKIKNMNST</p> <p>LTFVTLGEL RARRAEDGIV LDLPLYPAHP QDFHEVEDLI KTAIGNTLVQ DICYSPDTQK</p> <p>LLVRLSDVYN RSFLENLKVN TENLLQVENT GKVKGLLITL KGEPGGQTQA FDFYSRYFAP</p> <p>WVGVAEDPVT GSAHAVLSSY WSQHLGKKEM HAFQCShRGG ELGISLRPDG</p> <p>RVDIRGGA AV VLEGLTA</p>		
Quality Control	Endotoxin: Less than 0.1 ng/µg (1 IEU/µg) as determined by LAL test.		
Directions for Use	Always centrifuge tubes before opening. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com

Product Data Sheet

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com