

## Recombinant Human PCSK9 Protein

Catalog #	Source	Reactivity	Applications
CRP2064	Human cells	Human	E, WB, SDS-PAGE, MS

**Description** Recombinant Human PCSK9 Protein is produced by our mammalian expression system and the target gene encoding Gln31-Gln692 is expressed with a 6His/HA/AVI tag at the C-terminus.

**Form** Liquid in a 0.2 μM filtered solution of 50mM HEPES, 150mM NaCl, 20% Glycerol, pH7.4.

**Gene Symbol** PCSK9

**Alternative Names** NARC1; Proprotein convertase subtilisin/kexin type 9; Neural apoptosis-regulated convertase 1; NARC-1; Proprotein convertase 9; PC9; Subtilisin/kexin-like protease PC9

**Entrez Gene** 255738 (Human)

**SwissProt** Q8NBP7 (Human)

**Purity** Greater than 95% as determined by reducing SDS-PAGE.

**Chemical Structure** QEDEDGDYEE LVLALRSEED GLAEAPEHGT TATFHRC AKD PWRLPGTYVV VLKEETHLSQ SERTARRLQA QAARRGYLTK ILHVFHGLLP GFLVKMSGDL LELALKLPHV DYIEEDSSVF AQSIPWNLER ITPPRYRADE YQPPDGGSLV EVYLLDTSIQ SDHREIEGRV MVTD FENVPE EDGTRFHRQA SKCD SHGTHL AGVVSGRDAG VAKGASM RSL RVLNCQ GKGT VSGTLIGLEF IRKSQLVQPV GPLVLLPLA GGYSRVLNAA CQLR ARAGVV LVTAAGNFRD DACLYSPASA PEVITVGATN AQDQPVT LGT LGTNFGRCVD LFAPGEDIIG ASSDCSTCFV SQSGTSQAAA HVAGIAAMML SAPELT LAE LRQLIH FSA KDVINEAWFP EDQRVLT PNL VAALPPSTHG AGWQLFCRTV WSAHSGPTRM ATAIARCAPD EELLSCSSFS RSGKRRGERM EAQGGKLVCR AHNAFGGEGV YAIARCCLLP QANCSVHTAP PAEASMGTRV HCHQQGHVLT

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com

# Product Data Sheet

GCSSHWEVED LGTHKPPVLR PRGQPNQCVG HREASIHASC CHAPGLECKV KEHGIPAPQE  
QVTVACEEGW TLTGCSALPG TSHVLGAYAV DNTCVVRSRD VSTTGSTSEE AVTAVAICCR  
SRHLAQASQE LQEFGGHHHH HHHHGYPYDV PDYASGGGLN DIFEAQKIEW HE

**Quality Control**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Directions for Use**

Always centrifuge tubes before opening. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

**Storage/Stability**

Store it at -20°C to -80°C for one year.

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

**COHESION BIOSCIENCES LIMITED**

**WEB**

[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**

[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)