

Recombinant Human PPCDC Protein

Catalog #	Source	Reactivity	Applications
CRP2133	E. coli	Human	E, WB, SDS-PAGE, MS
Description	Recombinant Human PPCDC Protein is produced by our E.coli expression system and the target gene encoding Met1-Ser204 is expressed with a 6His tag at the N-terminus.		
Form	Liquid in a 0.2 μ M filtered solution of 20mM Tris-HCl, 50mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0.		
Gene Symbol	PPCDC		
Alternative Names	COAC; Phosphopantothenoylcysteine decarboxylase; PPC-DC; CoaC		
Entrez Gene	60490 (Human)		
SwissProt	Q96CD2 (Human)		
Purity	Greater than 95% as determined by reducing SDS-PAGE.		
Chemical Structure	MGSSHHHHHH SSGLVPRGSH MEPKASCPAA APLMERKFHV LVGVTGSVAA LKLPLLVSKL LDIPGLEVAV VTTERAKHFY SPQDIPVTLY SDADEWEMWK SRSDPVLHID LRRWADLLL APLDANTLGK VASGICDNL TCMRAWDRS KPLLFPCAMN TAMWEHPITA QQVDQLKAFG YVEIPCAKK LVCGDEGLGA MAEVGTIVDK VKEVLFQHSG FQQS		
Quality Control	Endotoxin: Less than 0.1 ng/ μ g (1 IEU/ μ g) as determined by LAL test.		
Directions for Use	Always centrifuge tubes before opening. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com