

## Recombinant Human PPM1A Protein

Catalog #	Source	Reactivity	Applications
CRP2139	E. coli	Human	E, WB, SDS-PAGE, MS
<b>Description</b>	Recombinant Human PPM1A Protein is produced by our E. coli expression system and the target gene encoding Gly2-Trp382 is expressed with a 6His tag at the C-terminus.		
<b>Form</b>	Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.		
<b>Gene Symbol</b>	PPM1A		
<b>Alternative Names</b>	PPPM1A; Protein phosphatase 1A; Protein phosphatase 2C isoform alpha; PP2C-alpha; Protein phosphatase 1A		
<b>Entrez Gene</b>	5494 (Human)		
<b>SwissProt</b>	P35813 (Human)		
<b>Purity</b>	Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.		
<b>Chemical Structure</b>	<p>MGAFLDKPKM EKHNAQGQGN GLRYGLSSMQ GWRVEMEDAH TAVIGLPSGL</p> <p>ESWSFFAVYD GHAGSQVAKY CCEHLLDHIT NNQDFKGSAG APSVENVKNG IRTGFLEIDE</p> <p>HMRVMSEKKH GADRSSTAV GVLISQHTY FINCGDSRGL LCRNRKVHFF TQDHPKSNPL</p> <p>EKERIQNAGG SVMIQRVNGS LAVSRALGDF DYKCVHGKGP TEQLVSPEPE VHDIERSEED</p> <p>DQFIIACDG IWDVMGNEEL CDFVRSRLEV TDDLEKVCNE VVDTCLYKGS RDNMSVILIC</p> <p>FPNAPKVSPE AVKKEAELDK YLECRVEEII KKQGEGVPDL VHVMRTLASE NIPSLPPGGE</p> <p>LASKRNVIEA VYNRLNPYKN DDTDSTSTDD MWLEHHHHHH</p>		
<b>Quality Control</b>	Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.		
<b>Directions for Use</b>	Always centrifuge tubes before opening. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		
<b>Storage/Stability</b>	Store it at -20°C to -80°C for one year.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com

# Product Data Sheet

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

**WEB**

[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**

[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)