

## Recombinant Human SLAMF8 Protein

Catalog #	Source	Reactivity	Applications
CRP2315	Human cells	Human	E, WB, SDS-PAGE, MS
<b>Description</b>	Recombinant Human SLAMF8 Protein is produced by our mammalian expression system and the target gene encoding Ala23-Asp233 is expressed with a 6His tag at the C-terminus.		
<b>Form</b>	Lyophilized from a 0.2 $\mu$ M filtered solution of PBS, pH 7.4.		
<b>Gene Symbol</b>	SLAMF8		
<b>Alternative Names</b>	BLAME; SLAM family member 8; B-lymphocyte activator macrophage expressed; BCM-like membrane protein; CD353		
<b>Entrez Gene</b>	56833 (Human)		
<b>SwissProt</b>	Q9P0V8 (Human)		
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE.		
<b>Chemical Structure</b>	<p>AQVLSKVGGS VLLVAARPPG FQVREAIWRS LWPSEELLAT FFRGSLETLY HSRFLGRAQL</p> <p>HSNLSLELGP LESGDSGNFS VLMVDTRGQP WTQTLQLKVY DAVPRPVVQV FIAVERDAQP</p> <p>SKTCQVFLSC WAPNISEITY SWRRETTMDF GMEPHSLFTD GQVLSISLGP GDRDVAYSCI</p> <p>VSNPVSWDLA TVTPWDSCHH EAAPGKASYK DVDHHHHHHH</p>		
<b>Quality Control</b>	Endotoxin: Less than 0.1 ng/ $\mu$ g (1 IEU/ $\mu$ g) as determined by LAL test.		
<b>Directions for Use</b>	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 $\mu$ g/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
<b>Storage/Stability</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)

# Product Data Sheet

Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

**WEB**

[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**

[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)