

Recombinant Human SULT1A2 Protein

Catalog #	Source	Reactivity	Applications
CRP2351	E. coli	Human	E, WB, SDS-PAGE, MS
Description	Recombinant Human SULT1A2 Protein is produced by our E. coli expression system and the target gene encoding Met1-Leu295 is expressed with a 6His tag at the N-terminus.		
Form	Liquid in a 0.2 μM filtered solution of 20mM Tris, 100mM NaCl, pH 8.0.		
Gene Symbol	SULT1A2		
Alternative Names	STP2; Sulfotransferase 1A2; ST1A2; Aryl sulfotransferase 2; Phenol sulfotransferase 2; Phenol-sulfating phenol sulfotransferase 2; P-PST 2		
Entrez Gene	6799 (Human)		
SwissProt	P50226 (Human)		
Purity	Greater than 95% as determined by reducing SDS-PAGE.		
Chemical Structure	<p>MGSSHHHHH SGLVPRGSH MELIQDISRP PLEYVKGVPL IKYFAEALGP LQSFQARPDD</p> <p>LLISTYPKSG TTWVSQILDM IYQGGDLEKC HRAPIFMRVP FLEFKVPGIP SGMETLKNTP</p> <p>APRLLKTHLP LALLPQTLTD QKVKVVYVAR NAKDVAVSYY HFYHMAKVYP HPGTWESFLE</p> <p>KFMAGEVSYG SWYQHVQEWV ELSRTHPVLY LFYEDMKENP KREIQKILEF VGRSLPEETV</p> <p>DLMVEHTSFK EMKKNPMTNY TTVRREFMDH SISPFMRKGM AGDWKTTFTV</p> <p>AQNERFDADY AEKMAGCSLS FRSEL</p>		
Quality Control	Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.		
Directions for Use	Always centrifuge tubes before opening. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com