

## Recombinant Mouse CD244 Protein

Catalog #	Source	Reactivity	Applications
CRP2892	Human cells	Mouse	E, WB, SDS-PAGE, MS
<b>Description</b>	Recombinant Mouse CD244 Protein is produced by mammalian expression system and the target gene encoding Gln20-Asn221 is expressed with a 6His tag at the C-terminus.		
<b>Form</b>	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.		
<b>Gene Symbol</b>	CD244		
<b>Alternative Names</b>	2B4; Natural killer cell receptor 2B4; NK cell activation-inducing ligand; NAIL; NK cell type I receptor protein 2B4; NKR2B4; h2B4; SLAM family member 4; SLAMF4; Signaling lymphocytic activation molecule 4; CD244		
<b>Entrez Gene</b>	18106 (Mouse)		
<b>SwissProt</b>	Q07763 (Mouse)		
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE.		
<b>Chemical Structure</b>	<p>QDCPDSSEEV VGVSGKPVQL RPSNIQTKDV SVQWKKTEQG SHRKIEILNW YNDGPSWSNV</p> <p>SFSDIYGFDY GDFALSIKSA KLQDSGHYLL EITNTGGKVC NKNFQLLLLD HVETPNLKAQ</p> <p>WKPWTNGTCQ LFLSCLVTKD DNVSALYRG STLISNQRNS THWENQIDAS SLHTYTCNVS</p> <p>NRASWANHTL NPTHGCSVP SNHHHHHH</p>		
<b>Quality Control</b>	Endotoxin: Less than 0.1 ng/µg (1 IEU/µg) as determined by LAL test.		
<b>Directions for Use</b>	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in ddH <sub>2</sub> O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
<b>Storage/Stability</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature		
<b>Application key:</b>	E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference		
<b>Species reactivity key:</b>	H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish		

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com

## Product Data Sheet

for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days.

Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**

[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**

[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)