

Recombinant Mouse CD307a Protein

| Catalog # | Source | Reactivity | Applications |
|--------------------------------|--|------------|---------------------|
| CRP2905 | Human cells | Mouse | E, WB, SDS-PAGE, MS |
| Description | Recombinant Mouse CD307a Protein is produced by mammalian expression system and the target gene encoding Ala17-Ser219 is expressed with a 6His tag at the C-terminus. | | |
| Form | Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4. | | |
| Gene Symbol | FCRL1 | | |
| Alternative Names | FCRH1; IFGP1; IRTA5; Fc receptor-like protein 1; FcR-like protein 1; FcRL1; Fc receptor homolog 1; FcRH1; IFGP family protein 1; hIFGP1; Immune receptor translocation-associated protein 5; CD307a | | |
| Entrez Gene | 229499 (Mouse) | | |
| SwissProt | Q8R4Y0 (Mouse) | | |
| Purity | Greater than 95% as determined by reducing SDS-PAGE. | | |
| Chemical Structure | AGISDVSLKT RPPGGWVMEG DKLVLICSDV RVTGNITYFW YRGALGFQLE TKTQPSTLAE FEISDMKQSD ADQYYCAAND GHDPIASELV SIHVRVPVSR PVLTFGDSGT QAVLGDLEVEL HCKALRGSPF IFYQFYHESI ILGNSSAPSG GGASFNFSLT AEHSGNFSCE ASNGQGAQRS EVVALNLTGL SLVPTENGIS HLSHHHHHH | | |
| Quality Control | Endotoxin: Less than 0.1 ng/µg (1 IEU/µg) as determined by LAL test. | | |
| Directions for Use | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in ddH ₂ O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. | | |
| Storage/Stability | Lyophilized protein should be stored at -20°C, though stable at room temperature | | |
| Application key: | E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference | | |
| Species reactivity key: | H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish | | |

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com

Product Data Sheet

for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days.

Aliquots of reconstituted samples are stable at -20°C for 3 months.

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com