

Recombinant Rat CD86 Protein

| Catalog # | Source | Reactivity | Applications |
|---------------------------|--|------------|---------------------|
| CRP3084 | Human cells | Rat | E, WB, SDS-PAGE, MS |
| Description | Recombinant Rat CD86 Protein is produced by mammalian expression system and the target gene encoding V1a29-Lys247 is expressed with a 6His tag at the C-terminus. | | |
| Form | Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4. | | |
| Gene Symbol | CD86 | | |
| Alternative Names | CD28LG2; T-lymphocyte activation antigen CD86; Activation B7-2 antigen; B70; BU63; CTLA-4 counter-receptor B7.2; FUN-1; CD86 | | |
| Entrez Gene | 56822 (Rat) | | |
| SwissProt | O35531 (Rat) | | |
| Purity | Greater than 95% as determined by reducing SDS-PAGE. | | |
| Chemical Structure | VPVKRQAYFN STAYLPCPFT KAQNISPSEL VVFWQDRKKS VLYEHYLGAE KLDNVNAKYL GRTSFDRDNQ ALRLHNVQIK DTGLYDCFIQ QKTPTGSIIIL QQWETELSVI ANFSEPEIEE AQNETRNTGI NLTCSSKQGY PKPTKMYFLI TNSTNEYGDN MQISQDNVTK LFSVSISLSL PFPDGVYNMT IVCILETESM NISSKPHNMV FSQPQFDRKH HHHHH | | |
| Quality Control | Endotoxin: Less than 0.1 ng/µg (1 IEU/µg) as determined by LAL test. | | |
| Directions for Use | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in ddH ₂ O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. | | |
| Storage/Stability | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. | | |

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com

Product Data Sheet

Aliquots of reconstituted samples are stable at -20°C for 3 months.

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com