

## Recombinant Rat IL-2 Protein

Catalog #	Source	Reactivity	Applications
CRP3090	Human cells	Rat	E, WB, SDS-PAGE, MS
<b>Description</b>	Recombinant Rat IL-2 Protein is produced by mammalian expression system and the target gene encoding Ala21-Gln155 is expressed with a 6His tag at the C-terminus.		
<b>Form</b>	Liquid in a 0.2 µm filtered solution of 5mM NaH <sub>2</sub> PO <sub>4</sub> , 5mM Citric acid, 150mM NaCl, pH4.0.		
<b>Gene Symbol</b>	IL2		
<b>Alternative Names</b>	Interleukin-2; IL-2; T-cell growth factor; TCGF; Aldesleukin		
<b>Entrez Gene</b>	116562 (Rat)		
<b>SwissProt</b>	P17108 (Rat)		
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE.		
<b>Chemical Structure</b>	APTSSPAKET QQHLEQLLLD LQVLLRGIDN YKNLKLPMML TFKFYLPKQA TELKHLQCLE NELGALQRVL DLTQSKSFHL EDAGNFISNI RVTVVVKLKGS ENKFECQFDD EPATVVEFLR RWIAICQSII STMTQH HHHHH H		
<b>Quality Control</b>	Endotoxin: Less than 0.1 ng/µg (1 IEU/µg) as determined by LAL test.		
<b>Directions for Use</b>	Please minimize freeze-thaw cycles.		
<b>Storage/Stability</b>	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**

[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**

[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**

[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**

[custom@cohesionbio.com](mailto:custom@cohesionbio.com)