


Anti-GFP-tag Antibody-AcalephFluor555 labeled

Catalog #	Source	Reactivity	Applications
CSA5405	Mouse	N/A	WB, IF, FC
Description	Mouse monoclonal antibody to GFP-tag-AcalephFluor555 labeled		
Immunogen	Recombinant protein corresponding to GFP-tag.		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes GFP and GFP tag fusion proteins.		
Clonality	Monoclonal		
Conjugation	AcalephFluor555		
Form	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide		
Dilution	WB (1/2000 - 1/5000), IF (1/500 - 1/2000), FC (1/1000 - 1/4000)		
Fluorophore	Ex/Em = 555/565 nm		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.		



Line colors represent the approximate visible colors of the wavelength maxima.

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com