

L-Lactate Dehydrogenase

Catalog #	Source	Reactivity	Applications
CZM1005	Porcine heart		
Description	L-Lactate Dehydrogenase from Porcine heart		
Specificity	>90 U/mg powder		
Form	White Lyophilized Powder		
Alternative Names	EC 1.1.1.27; L-lactic acid dehydrogenase; L-lactic dehydrogenase		
CAS Number	9001-60-9		
Molecular Weight	35 kD		
Purity	>90%		
Storage/Stability	Shipped at 4°C. Store at -20°C for 2 years.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com